

Product catalogue
HAEMODIALYSIS

Your guide to Vantive
Haemodialysis products

Australia and New Zealand





TABLE OF CONTENTS

Haemodialysis Systems	4
Disinfection	7
Dialysers	8
ArtiSet Bloodlines	13
Novaline Bloodlines	16
Concentrates	21
Needles	26
Dressings	29
Water Systems	30


HAEMODIALYSIS SYSTEMS

Vantive HD portfolio with the **AK 98** system is designed to address key healthcare challenges

The **AK 98** system is purposefully designed to **simplify workflows** and **deliver operational efficiency**

Suitable for the majority of HD patients, including those requiring removal of large-middle molecules

Small and compact design contributes to **reducing the environment impact**



**Lower HD
Program Costs**



**Improve
Sustainability**



**Enable Simple
Usability**



**Deliver Clinical
Outcomes**

HAEMODIALYSIS SYSTEMS

AK 98 V3 BIO

Fixed screen solution perfect for satellite and regional units. Delivers **HDx** and HD therapies in both double and single-needle modes.



Product Code: 955603

Compatibility: **SoftPac** Citrate + **BiCart** Concentrate System
CleanCart A & C for Disinfection

Pack Factor: 1

Standard Features:

- Graphical User Interface
- **Diascan** (Real-Time Kt/V)
- Black Box Data Recording
- Blood Pressure Monitoring (BPM)
- Auto Prime
- Auto Rinseback
- UF and Sodium Profiling
- Concentrate Standby
- Water Saving Mode
- Battery Backup
- IT Connectivity via HL7 Messaging

AK 98 V3 SELF CARE

Flexible screen positioning is perfect for self-care and home dialysis patients. Delivers **HDx** and HD therapies in both double and single-needle modes.



Product Code: 955604

Compatibility: **SoftPac** Citrate + **BiCart** Concentrate System
CleanCart A & C for Disinfection

Pack Factor: 1

Standard Features:

- Graphical User Interface
- **Diascan** (Real-Time Kt/V)
- Black Box Data Recording
- Blood Pressure Monitoring (BPM)
- Auto Prime
- Auto Rinseback
- UF and Sodium Profiling
- Concentrate Standby
- Water Saving Mode
- Battery Backup
- IT Connectivity via HL7 Messaging

ULTRAFILTERS

For the purification of incoming water for dialysis fluid as well as purification of dialysis fluid to obtain microbiological high quality dialysis fluid, minimising the risk of exposure to bacteria and endotoxins.



U9000 PLUS


Product Code: 955825
Compatibility: Artis Physio Plus
Artis Physio
AK 98 / AK 96
Shelf Life: 24 Months
Pack Factor: 16

U9000 PLUS FOR HOME PATIENTS

Product Code: 955825X1
Compatibility: AK 98 / AK 96
Shelf Life: 24 Months
Pack Factor: 1

DISINFECTION

CLEANCART CARTRIDGES



CLEANCART C

For the preparation of a citric acid solution used for removing calcium and magnesium from the dialysis machine's fluid circuit. **CleanCart C** must be used in combination with the heat disinfection program of the corresponding dialysis machine.

Contains 32g of citric acid anhydrous powder.

Product Code:	955851
Compatibility:	Artis Physio Plus / Artis Physio AK 98 / AK 96 / AK 95 S AK 200 Ultra S / AK 200 S
Shelf Life:	24 Months
Pack Factor:	33



CLEANCART A

For the preparation of a sodium carbonate solution used for removing organic deposits, fats and proteins from the dialysis machine's fluid circuit. **CleanCart A** must be used in combination with the heat disinfection program of the corresponding dialysis machine.

Contains 13g of anhydrous sodium carbonate powder.

Product Code:	955850
Compatibility:	Artis Physio Plus / Artis Physio AK 98 / AK 96 / AK 95 S AK 200 Ultra S / AK 200 S
Shelf Life:	24 Months
Pack Factor:	33

DIALYSERS

POLYFLUX L

Low flux dialyser featuring a **Polyamix** membrane.



POLYFLUX 14L

Product Code:	102057
Compatibility:	HD
Surface Area:	1.4m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	500ml

POLYFLUX 17L

Product Code:	102058
Compatibility:	HD
Surface Area:	1.7m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	500ml

POLYFLUX 21L

Product Code:	104176
Compatibility:	HD
Surface Area:	2.1m ²
Shelf Life:	24 Months
Pack Factor:	24

POLYFLUX H

High flux dialyser featuring a **Polyamix** membrane.



POLYFLUX 140H

Product Code:	103530
Compatibility:	HD / HDF
Surface Area:	1.4m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	500ml



POLYFLUX 170H

Product Code:	103579
Compatibility:	HD / HDF
Surface Area:	1.7m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	500ml



POLYFLUX 210H

Product Code:	115821
Compatibility:	HD / HDF
Surface Area:	2.1m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	500ml

POLYFLUX 2H/6H

High flux dialyser featuring a **Polyamix** membrane for low body weight patients, typically children.



POLYFLUX 2H

Product Code:	112468
Compatibility:	HD / HDF
Surface Area:	0.2m ²
Shelf Life:	24 Months
Pack Factor:	16
Min Recommended Priming Volume:	500ml



POLYFLUX 6H

Product Code:	103403
Compatibility:	HD / HDF
Surface Area:	0.6m ²
Shelf Life:	24 Months
Pack Factor:	16
Min Recommended Priming Volume:	1,000ml

THERANOVA

A unique **MCO** dialyser that efficiently targets the removal of large-middle molecules, enabling **HDx** therapy.



THERANOVA 400

Product Code:	115578
Compatibility:	HDx
Surface Area:	1.7m ²
Shelf Life:	24 Months
Pack Factor:	24

THERANOVA 500

Product Code:	115579
Compatibility:	HDx
Surface Area:	2.0m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	300ml

THERALITE

An **HCO** dialyser designed to target the removal of free light chain (FLC) proteins and other plasma components with molecular weights up to 45 kDa.



THERALITE

Product Code:	113092
Compatibility:	HD
Surface Area:	2.1m ²
Shelf Life:	24 Months
Pack Factor:	1
Min Recommended Priming Volume:	1,000ml

EVODIAL

A specialised dialyser for patients with a high risk of bleeding. **HeprAN** (heparin-grafted **AN69 ST**) membrane for patients requiring reduced or heparin-free dialysis.



EVODIAL 1.0

Product Code:	110654A
Compatibility:	HD
Surface Area:	1.05m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	1,000ml (at UFR = 2,000ml / h)

EVODIAL 1.3

Product Code:	110653A
Compatibility:	HD
Surface Area:	1.30m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	1,000ml (at UFR = 2,000ml / h)

EVODIAL 1.6

Product Code:	110652A
Compatibility:	HD
Surface Area:	1.65m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	1,000ml (at UFR = 2,000ml / h)

EVODIAL 2.2

Product Code:	110651A
Compatibility:	HD
Surface Area:	2.15m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	1,000ml (at UFR = 2,000ml / h)

NEPHRAL ST

A dialyser specialised for high biocompatibility and adsorption - featuring a heparin adsorptive hydrogel (**AN69 ST**) membrane for reduced use of heparin during dialysis and designed to provide effective removal of uremic toxins and inflammatory mediators by adsorption. **Nephral ST** is a dialyser for acetate-free biofiltration and may be an alternative for patients experiencing hypersensitivity reactions to conventional membranes, and may also help enhance removal of cytokines such as TNF- α , IL-6 and IL-8.




NEPHRAL ST 200

Product Code:	101939A
Compatibility:	HD
Surface Area:	1.05m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	1,000ml (at UFR = 2,000ml / h)
Heparinised Solution:	5,000 IU / L

NEPHRAL ST 500

Product Code:	101942A
Compatibility:	HD
Surface Area:	2.15m ²
Shelf Life:	24 Months
Pack Factor:	24
Min Recommended Priming Volume:	1,000ml (at UFR = 2,000ml / h)
Heparinised Solution:	5,000 IU / L


ARTISET BLOODLINES



ARTISET HD DNL HC

Integrated blood circuit for double needle (DNL) haemodialysis.


Product Code:	955075
Compatibility:	Artis Physio Plus Artis Physio
Sterilisation:	Gamma / Beta Radiation
Shelf Life:	36 Months
Pack Factor:	20



ARTISET SN HC

Integrated blood circuit for single needle (SN) haemodialysis.

Product Code:	114533
Compatibility:	Artis Physio Plus Artis Physio
Sterilisation:	Gamma Radiation
Shelf Life:	36 Months
Pack Factor:	7




ARTISET PREPOST

An integrated blood circuit for pre and post dilution haemodialysis.

Product Code:	955077
Compatibility:	Artis Physio Plus Artis Physio
Sterilisation:	Gamma Radiation
Shelf Life:	36 Months
Pack Factor:	20

ARTISET BLOODLINE ACCESSORIES

All accessories are Latex free.




DEHP-FREE

ULTRALINE HD

Provides the capability to use online fluids for priming and substitution fluid. Used in conjunction with **ArtiSet** DNL or SN bloodlines on the **Artis Physio Plus** system.

Product Code:	955600
Compatibility:	Artis Physio Plus
Sterilisation:	Beta Radiation
Shelf Life:	36 Months
Pack Factor:	80




DEHP-FREE

ULTRALINE HDF

Provides the capability to use online fluids for priming and substitution fluid. Used in conjunction with **ArtiSet** bloodlines on the **Artis Physio Plus** system.

Product Code:	955599
Compatibility:	Artis Physio Plus
Sterilisation:	Beta Radiation
Shelf Life:	36 Months
Pack Factor:	40




DEHP-FREE

ULTRA PRIME LINE

Provides the capability to use online fluids for priming - fitted with a sterilising membrane and a T-connector for rinseback. Used in conjunction with **ArtiSet** bloodlines on the **Artis Physio** system.

Product Code:	955490
Compatibility:	Artis Physio
Sterilisation:	Beta Radiation
Pack Factor:	56




DEHP-FREE

ULTRA HDF POST LINE

A single-use ultrafilter dedicated for post-dilution modalities. Used in conjunction with the **ArtiSet** Prepost bloodline on the **Artis Physio** system.

Product Code:	955037
Compatibility:	Artis Physio
Sterilisation:	Beta Radiation
Shelf Life:	12 Months
Pack Factor:	20

ARTISET BLOODLINE ACCESSORIES CONT'D




DEHP-FREE

ULTRA HDF LINE

A single-use ultrafilter dedicated for pre-dilution modalities. Used in conjunction with the **ArtiSet** Prepost or **ArtiSet** HD DNL HC bloodlines on the **Artis Physio** system.

Product Code:	115283
Compatibility:	Artis Physio
Sterilisation:	Gamma Radiation
Shelf Life:	12 Months
Pack Factor:	18



DEHP-FREE

PRIMING ADAPTOR

Used as a recirculation adaptor during treatment to connect the arterial and venous patient lines or to connect the blood tubing set to the waste handling option (WHO) during priming.

Product Code:	106283
Compatibility:	Artis Physio Plus / Artis Physio AK 98 / AK 96 / AK 95 S AK 200 Ultra S / AK 200 S
Sterilisation:	Gamma
Shelf Life:	36 Months
Pack Factor:	100



DEHP-FREE


GMB-SP06

500ml sampling bag for micro-biological or chemical dialysis fluid sampling. Male connection.

Product Code:	115648
Compatibility:	Artis Physio Plus / Artis Physio AK 200 S / AK 200 ULTRA S
Sterilisation:	EtO
Shelf Life:	36 Months
Pack Factor:	36

NOVALINE BLOODLINES

NOVALINE BLOODLINES *FOR HD ON AK SYSTEMS*




BL 14

Double Needle Treatment

Product Code:	955304
Compatibility:	AK 98 / AK 96 / AK 95 S / AK 200 S
Sterilisation:	Radiation
DEHP-Free:	Yes
Shelf Life:	36 Months
Pack Factor:	24
Blood Pathway Volume:	133ml*
Pump Segment:	8 x 12mm - 250mm
Infusion Line:	Yes
Pressure Transducer:	1
Pre-Dialyser Chamber:	No
Dialyser Connector:	Dual
Removable Labelling:	No

NOVALINE BLOODLINES *FOR OTHER HD SYSTEMS*




BL 37

Double Needle Treatment

Product Code:	955459
Compatibility:	Fresenius 4008
Sterilisation:	Radiation
DEHP-Free:	Yes
Shelf Life:	36 Months
Pack Factor:	20
Blood Pathway Volume:	181ml*
Pump Segment:	8 x 12mm - 355mm
Infusion Line:	Yes
Pressure Transducer:	1
Pre-Dialyser Chamber:	Yes
Dialyser Connector:	Dual
Removable Labelling:	No
BVS + Combi Connector:	No

* +/- 10%

NOVALINE BLOODLINES FOR PAEDIATRIC AND LOW / VERY LOW VOLUME PATIENTS ON AK SYSTEMS



* +/- 10%

BL 122P

Double Needle Treatment

Product Code: 955457

Compatibility: AK 98 / AK 96 / AK 95 S / AK 200
ULTRA S / AK 200 S

Sterilisation: Radiation

DEHP-Free: Yes

Shelf Life: 36 Months

Pack Factor: 24

Blood Pathway Volume: 116ml*

Pump Segment: 6.35 x 10.35mm - 250mm

Infusion Line: No

Pressure Transducer: 1


Pre-Dialyser Chamber: Yes

Dialyser Connector: Dual

Removable Labelling: Yes

BVS + Combi Connector: No


NOVALINE BLOODLINES FOR PEDIATRIC AND LOW / VERY LOW VOLUME PATIENTS ON AK SYSTEMS CONT'D



BL 121P HDF

Double Needle Treatment

Product Code:	955464	Blood Pathway Volume:	87ml*
Compatibility:	AK 98 / AK 200 S AK 200 ULTRA S	Pump Segment:	6.35 x 10.35mm - 247mm
Sterilisation:	Radiation	Infusion Line:	No
DEHP-Free:	Yes	Pressure Transducer:	1
Shelf Life:	36 Months	Pre-Dialyser Chamber:	No
Pack Factor:	24	Dialyser Connector:	Dual
		Removable Labelling:	Yes
		BVS + Combi Connector:	Yes



BL 120N


Double Needle Treatment

Product Code:	955463	Blood Pathway Volume:	32ml*
Compatibility:	AK 98 / AK 200 S AK 200 ULTRA S	Pump Segment:	3.9 x 7.9mm - 250mm
Sterilisation:	Radiation	Infusion Line:	No
DEHP-Free:	Yes	Pressure Transducer:	1
Shelf Life:	36 Months	Pre-Dialyser Chamber:	No
Pack Factor:	36	Dialyser Connector:	Dual
		Removable Labelling:	Yes
		BVS + Combi Connector:	No

* +/- 10%

NOVALINE BLOODLINE ACCESSORIES


All accessories are Latex free.



C705

Expansion chamber for emergency switch from double needle to single needle / single pump mode.


Product Code:	101739
Compatibility:	AK 98 / AK 200 S AK 200 ULTRA S
Sterilisation:	EtO
Shelf Life:	24 Months
Pack Factor:	32



PRIMING ADAPTOR

Used as a recirculation adaptor during treatment to connect the arterial and venous patient lines.


Product Code:	106283
Compatibility:	AK 98 / AK 96 / AK 95 S AK 200 S / AK 200 ULTRA S
Sterilisation:	Gamma
Shelf Life:	36 Months
Pack Factor:	100



S-660-C

Y connector (10cm lines with male / 2 x female connections) to use with single needle dialysis with a single lumen catheter.

Product Code:	101354
Compatibility:	All
Sterilisation:	EtO
Shelf Life:	36 Months
Pack Factor:	50




S-1016

Extension line (20cm) with male / female luer to extend the length of the patient line.

Product Code:	101023
Compatibility:	All
Sterilisation:	EtO
Shelf Life:	36 Months
Pack Factor:	50


NOVALINE BLOODLINE ACCESSORIES *CONT'D*



SP-C35

Pressure transducer protector (rigid connection, with male / female connection) for use when only the protector needs to be changed or added.

Product Code:	106024
Compatibility:	AK 98 / AK 96 / AK 95 S AK 200 S / AK 200 ULTRA S
Sterilisation:	EtO
Shelf Life:	36 Months
Pack Factor:	70



SP-C13

Unvented spike with female connection for use with saline bags.

Product Code:	101593
Compatibility:	All
Sterilisation:	EtO
Shelf Life:	24 Months
Pack Factor:	75

SELECTBAG CITRATE



Pack Factor: 6

21

SOFTPAC CITRATE



SOFTPAC CITRATE

An acetate-free concentrate with all the beneficial properties of citrate, to be used with a dry sodium bicarbonate cartridge (**BiCart**) for online preparation of haemodialysis fluid.

Dilution ration 1+44

Compatibility: AK 98 / AK 96 / AK 95 S
AK 200 Ultra S / AK 200 S
Artis Physio / Artis Physio Plus

Shelf Life: 24 Months

Pack Factor 3.5L: 3

Pack Factor 5L: 2

PRODUCT CODE	PRODUCT NAME	Na ⁺ MMOL/L	HCO ₃ ⁻ MMOL/L	K ⁺ MMOL/L	Ca ²⁺ MMOL/L	Mg ²⁺ MMOL/L	Cl ⁻ MMOL/L	Glucose G/L
115753	SoftPac Citrate C140 3.5L	140	34	1	1.41	0.5	107.8	1
115755	SoftPac Citrate C165 3.5L	140	34	1	1.66	0.5	108.3	1
115760	SoftPac Citrate C240 3.5L	140	34	2	1.41	0.5	108.8	1
115761	SoftPac Citrate C250 3.5L	140	34	2	1.51	0.5	109.0	1
115762	SoftPac Citrate C265 3.5L	140	34	2	1.66	0.5	109.3	1
115763	SoftPac Citrate C275 3.5L	140	34	2	1.76	0.5	109.5	1
115767	SoftPac Citrate C340 3.5L	140	34	3	1.41	0.5	109.8	1
115768	SoftPac Citrate C350 3.5L	140	34	3	1.51	0.5	110.0	1
115769	SoftPac Citrate C365 3.5L	140	34	3	1.66	0.5	110.3	1
115770	SoftPac Citrate C375 3.5L	140	34	3	1.76	0.5	110.5	1
115774	SoftPac Citrate C450 3.5L	140	34	4	1.51	0.5	111.0	1
115775	SoftPac Citrate C465 3.5L	140	34	4	1.66	0.5	111.3	1
115824	SoftPac Citrate G150 5L	140	34	1	1.51	0.5	108.0	1
115825	SoftPac Citrate G165 5L	140	34	1	1.66	0.5	108.3	1
115784	SoftPac Citrate G250 5L	140	34	2	1.51	0.5	109.0	1
115785	SoftPac Citrate G265 5L	140	34	2	1.66	0.5	109.3	1
115786	SoftPac Citrate G275 5L	140	34	2	1.76	0.5	109.5	1
115789	SoftPac Citrate G350 5L	140	34	3	1.51	0.5	110.0	1
115790	SoftPac Citrate G365 5L	140	34	3	1.66	0.5	110.3	1
115791	SoftPac Citrate G375 5L	140	34	3	1.76	0.5	110.5	1

SOFTPAC



SOFTPAC

An acid concentrate for bicarbonate based haemodialysis and intended to be used with a dry sodium bicarbonate cartridge (**BiCart**) for online preparation of haemodialysis fluid.

Dilution ration 1+44

Compatibility: AK 98 / AK 96 / AK 95 S
AK 200 Ultra S / AK 200 S
Artis Physio / Artis Physio Plus

Shelf Life: 24 Months

Pack Factor 3.5L: 3

Pack Factor 5L: 2

PRODUCT CODE	PRODUCT NAME	Na ⁺	HCO ₃ ⁻	K ⁺	Ca ²⁺	Mg ²⁺	Cl ⁻	ACETATE	GLUCOSE
		MMOL/L	MMOL/L	MMOL/L	MMOL/L	MMOL/L	MMOL/L	MMOL/L	G/L
114933	SoftPac C197 3.5L	140	34	1	1.25	0.5	107.5	3	1
114941	SoftPac C298 3.5L	140	34	2	1.25	0.5	108.5	3	1
114948	SoftPac C395 3.5L	140	34	3	1.25	0.5	109.5	3	1
114952	SoftPac G194 5L	140	34	1	1.50	0.5	108.0	3	1
114961	SoftPac G298 5L	140	34	2	1.25	0.5	108.5	3	1
114965	SoftPac G394 5L	140	34	3	1.50	0.5	110.0	3	1

BICART SELECT



BICART SELECT COMBI PACK

The **BiCart** and **SelectCart** products are intended for preparation of haemodialysis solutions. The **BiCart** and **SelectCart** should be used together with **SelectBag** Citrate.

BiCart cartridge - 720g of NaHCO₃

SelectCart cartridge - 1200g NaCl

Product Code: 955848

Compatibility: **Artis Physio Plus**
Artis Physio
AK 200 S, AK 200 S ULTRA

Shelf Life: 24 Months

Pack Factor: 4

BICART

A powdered bicarbonate concentrate to be used with **SoftPac** or **SoftPac** Citrate for the preparation of online haemodialysis fluid.



BICART 720G

Product Code: 955834

Compatibility: **AK 98 / AK 96 / AK 95 S**
AK 200 S / AK 200 ULTRA S

(**Artis Physio Plus** and **Artis Physio** when used with **SoftPac** or **SoftPac** Citrate)

Shelf Life: 24 Months

Pack Factor: 10



BICART 1150G

Product Code: 955840

Compatibility: **AK 98 / AK 96 / AK 95 S**
AK 200 S / AK 200 ULTRA S


(**Artis Physio Plus** and **Artis Physio** when used with **SoftPac** or **SoftPac** Citrate)

Shelf Life: 24 Months

Pack Factor: 6

NEEDLES

AV FISTULA NEEDLES




AV FISTULA NEEDLES

Shelf Life: 24 Months

Pack Factor: 50 / Box

PRODUCT CODE	NEEDLE OD (G)	NEEDLE LENGTH (MM)	TUBE LENGTH (MM)	WING FIXED / ROTATING	SHARP / DULL (BUTTONHOLE)	STERILISATION
F14BSG	14	25	300	Rotating	Sharp	Gamma Radiation
F15BSG	15	25	300	Rotating	Sharp	Gamma Radiation
F16BSG	16	25	300	Rotating	Sharp	Gamma Radiation
F17BSG	17	25	300	Rotating	Sharp	Gamma Radiation
F15ASRG	15	20	150	Rotating	Sharp	Gamma Radiation
F15BSLG	15	32	300	Rotating	Sharp	Gamma Radiation

AV FISTULA NEEDLES - UNIPUNCTURE (SINGLE NEEDLE)



AV FISTULA NEEDLES - UNIPUNCTURE (SINGLE NEEDLE)

Shelf Life: 24 Months

Pack Factor: 40 / Box

PRODUCT CODE	NEEDLE OD (G)	NEEDLE LENGTH (MM)	TUBE LENGTH (MM)	WING FIXED / ROTATING	SHARP / DULL (BUTTONHOLE)	STERILISATION
FSN15ASG	15	25	150	Rotating	Sharp	Gamma Radiation
FSN16ASG	16	25	150	Rotating	Sharp	Gamma Radiation

PLUME SAFETY NEEDLE



PLUME SAFETY NEEDLE

Shelf Life: 24 Months

Pack Factor: 50 / Box

PRODUCT CODE	NEEDLE OD (G)	NEEDLE LENGTH (MM)	TUBE LENGTH (MM)	WING FIXED / ROTATING	SHARP / DULL (BUTTONHOLE)	STERILISATION
A15L30SGP	15	25	300	Rotating	Sharp	Electron Beam Radiation
A15M30SGP	15	20	300	Rotating	Sharp	Electron Beam Radiation
A16M30SGP	16	20	300	Rotating	Sharp	Electron Beam Radiation
A17M30SGP	17	20	300	Rotating	Sharp	Electron Beam Radiation
A16L30SGP	16	25	150	Rotating	Sharp	Electron Beam Radiation

THIRD PARTY AV NEEDLES

NIPRO NEEDLES

Shelf Life: 24 Months

PRODUCT CODE	NEEDLE OD (G)	NEEDLE LENGTH (MM)	TUBE LENGTH (MM)	WING FIXED / ROTATING	SHARP / DULL (BUTTONHOLE)
113144	14	25	400	Rotating	Sharp
114329	15	20	300	Rotating	Sharp
113145	15	25	400	Rotating	Sharp
113911	15	20	400	Rotating	Buttonhole
DAV525HTCW40N	15	25	400	Rotating	Buttonhole
2027717	15	20	400	Rotating	Sharp
109957	14	25	400	Rotating	Buttonhole

NEEDLE COVERS

NEEDLE COVERS

Plastic cover that converts standard AVF needles to a safety needle.

Product Code: 112075

Pack Factor: 1200 / Box

DRESSINGS



TIPSTOP

A compression dressing to stop bleeding after vascular punctures, it has a lyophilised alginate moulded into a polyamide compressive pad. It has a hypoallergenic adhesive film and a rigid PET protective housing.

Product Code:	956008
Sterilisation:	Gamma
Shelf Life:	24 Months
Pack Factor:	200 / Box



BAXTER PLASTER

Adhesive pressure bandage.

Product Code:	956009
Sterilisation:	EtO
Shelf Life:	24 Months
Pack Factor:	1,000 / Box

WATER SYSTEMS

CWP 800



WRO 300 H



CWP ONLINE

CWP 800

PRODUCT CODE	Product Name	Description	PACK FACTOR	Compatibility
115641	CWP 800 RO1 415V	Single RO device 1 x membrane, Pure water output 25l/min	1	
115642	CWP 800 RO2 415V	Single RO device 2 x membranes, Pure water output 48l/min	1	
115643	CWP 800 Ext 1	Extension device 1 x membrane, RO2 + Ext 1 Pure water output 66l/min	1	CWP 800 RO2 415V
115644	CWP 800 Ext 2	Extension device 1 x membrane, RO2 + Ext 2 Pure water output 78l/min	1	CWP 800 RO2 415V
115645	CWP 800 HW 415V	Hot water device effective tank volume 330L	1	CWP 800 RO1 CWP 800 RO2
115942	CWP 800 Clinic Panel	Remote user interface panel	1	CWP 800 RO1 CWP 800 RO2

WRO 300 H

Product Code:	955686	Product Description:
Pack Factor:	1	Single patient RO unit - Automated procedures minimise patient and caregiver involvement - Compatible dialysis machines can start and stop the WRO 300 H unit - Automatic shut-off at completion of disinfection process - The WRO 300 H unit can be programmed to automatically start heat disinfection at specified intervals - Provide 1.2l/min of permeate - Hot water disinfection - Low noise level

ACCESSORIES

PRODUCT CODE	Product Name	Description	PACK FACTOR	Compatibility
5167IN	HiSense Serim ULTRA 0.1 total chlorine test strips (1 x 100)	Total residual chlorine test strips (1 x 100)	1	
MANA91243	Aquadur Sensitive	Total water hardness test strips	1	
5105IN	Serim residual peroxide test strips	Residual peroxide test strips	1	Dialox 5L
DIALOX5KG	Dialox	Dialox 5L	1	

CWP ONLINE

The cloud-based, on-demand service for remote monitoring of your **Vantive** Water Systems

Comprehensive monitoring for informed decision making:

- Monitor your devices on the go
- Receive personalised event notifications
- Analyse data and get retrospective insights
- Achieve operational efficiency



Australia
Vantive Pty Ltd

1 Baxter Drive
Old Toongabbie
NSW 2146

 1800 920 133

 pro.vantive.com.au

New Zealand
Vantive Limited

33 Vestey Drive
Mt Wellington
Auckland 1060

 0800 436 770

 pro.vantive.com.au

Vantive, Gambro, AK 200 S, AK 200 ULTRA S, AK 95 S, AK 96, AK 98, AN69, ArtiSet, Artis Physio, Artis Physio Plus, BiCart, BiCart Select, CleanCart, CWP 800, Diascan, Evodial, HCO, HeprAN, HDx, MCO, Nephral, Novaline, Plume, Polyamix, Polyflux, SelectBag, SelectCart, SoftPac, Theranova, Theralite, TipStop, U9000 Plus and WRO 300 H are trademarks of Vantive Health LLC or its affiliates.

Baxter is a trademark of Baxter International Inc. or its subsidiaries.

Any other trademarks and product names appearing herein are the property of their respective owners.